


TheraTogs Orthotic Garment Systems

Sizing Chart

Sizing Guidelines



	Size Group	US – POUNDS / INCHES			METRIC – KG / CM			PTA System	Full Body System	Lower Extremity System
		Weight	STA	ATK	Weight	STA	ATK			
PEDIATRICS	Infant	15 – 25	8 – 10	7 - 10	7 - 11	20 – 25	18 - 25	☐ GPTA110	☐ GFB100	☐ GLE101
	Pre-School	25 – 45	11 - 13	11 - 12	12 - 20	28 - 33	28 - 31	☐ GPTA210	☐ GFB200	☐ GLE201
	Pediatric	45 – 65	14 - 15	13 - 18	21 - 30	36 – 38	33 - 46	☐ GPTA310	☐ GFB300	☐ GLE301
	Junior	65 – 110	16 - 17	19 - 21	31 - 54	41 - 44	48 - 53	☐ GPTA410	☐ GFB400	☐ GLE401
ADULT	Female Med	110 – 170	18 – 20	22 - 24	55 - 77	46 - 51	56 - 61	☐ GPTA610	☐ GFB600	☐ GLE601
	Female Lg	170 – 200	21 - 22	25 - 27	78 - 91	53 – 56	63 - 69	☐ GPTA620	☐ GFB620	☐ GLE621
	Male Med	110 – 170	18 – 19	22 - 24	55 - 77	46 - 49	56 - 61	☐ GPTA710	☐ GFB700	☐ GLE701
	Male Large	170 – 190	20 - 22	25 - 27	78 - 86	50 - 56	63 - 69	☐ GPTA720	☐ GFB720	☐ GLE721

	Size Name	Age (Mo)	US - Inches		Metric - CM		Wunzi Model #
			[A] Max Torso	[B] Max Belly	[A] Max Torso	[B] Max Belly	
WUNZI	Wunzi XX-Small	0 – 3	16	16	41	41	☐ GWZ25
	Wunzi X-Small	3 – 6	20	17	51	44	☐ GWZ50
	Wunzi Small	6 – 12	22	18	56	46	☐ GWZ75
	Wunzi Medium	12 – 18	24	20	61	51	☐ GWZ100
	Wunzi Large	18 – 24	27	21	69	53	☐ GWZ125
	Limb Kit, Small	0 – 12					☐ GWZLKS
	Limb Kit, Large	12 – 24					☐ GWZLKL

[A] Wunzi Torso Length



C-7 to navel over diaper

[B] Belly Girth




Around belly at navel

TherAlign Dragonfly TLSO


STEP 1: Record measurements

Garment Subsystem



Dragonfly Brace

Spine Length C-7 to Coccyx



STEP 2: Size & Select Garment Subsystem

SIZE GROUP	US – POUNDS / INCHES			METRIC – KG / CM			GARMENT SUBSYSTEM
	WEIGHT	STA	ATK	WEIGHT	STA	ATK	
Infant	15 – 25	8 – 10	7 - 10	7 - 11	20 – 25	18 - 25	☐ 110
Pre-School	26 – 45	11 - 13	11 - 12	12 - 20	28 - 33	28 - 31	☐ 210
Pediatric	46 – 65	14 - 15	13 - 18	21 - 30	36 – 38	33 - 46	☐ 310
Junior	66 – 120	16 - 17	19 - 21	31 - 54	41 - 44	48 - 53	☐ 410
Female Med.	121 – 170	18 – 20	22 - 24	55 - 77	46 - 51	56 - 61	☐ 610
Female Large	171 – 200	21 - 22	25 - 27	78 - 91	53 – 56	63 - 69	☐ 620
Male Med.	121 – 170	18 – 19	22 - 24	55 - 77	46 - 49	56 - 61	☐ 710
Male Large	171 – 190	20 - 22	25 - 27	78 - 86	50 - 56	63 - 69	☐ 720

STEP 3: Size & Select Dragonfly Brace

SPINE (IN)	SPINE (CM)	MODEL
10 - 12	26 - 32	☐ DFLY30
13 - 14	33 - 35	☐ DFLY33
15 - 16	40 - 43	☐ DFLY40
17 - 19	48 - 52	☐ DFLY48
20 - 21	53 - 56	☐ DFLY53

STEP 4: Measure Chest & Hip Girth

The Dragonfly normally ships with standard-length Rib and Pelvis straps.

Please measure and report Chest and Hip girth in the boxes at right.

We will substitute longer or shorter straps if indicated by your measurements.

Rib & Pelvis Strapping

